


Owner's Manual

Olympic Elite Exam Lift Table™

Model 51201, 51202

Serial # _____



9311 SE 36th Street, Unit 220
Mercer Island, WA 98040

Phone: 800-789-6867 • **Fax** 206-905-0839

Web: www.OlympicVeterinary.com

Email: Info@OlyVet.com

About this Manual

This manual contains information on the setup, use and upkeep of the Olympic Elite Exam Lift Table and should be read in its entirety before first use to assure safe operation and optimal performance.

This lift has a serial number you will need to reference when seeking technical support. This is stamped on the bottom right corner of the name plate, located on the side of the table, just under the Corian® top. For future reference, it is recommended you write this serial number on the front of this manual.

Table of Contents

Component Identification.....	3
Introduction	4
Safety	4
Unpacking and Setup	4
Assembly.....	5
Operation	5
Loading/Unloading Patients	5
Height Adjustment	5
Scale	5
Operating Procedure	5
Power	6
Moving the Elite Exam Lift Table.....	6
Cleaning	6
Routine Maintenance	6
Specifications.....	6
Warranty.....	7
Contact Us	7

List of Figures

Figure 1 – Unpacking the Exam Lift Table	3
Figure 2 – Component Identification.....	4
Figure 3 – Scale Control Panel.....	5

Component Identification

The following illustration demonstrates the key components of the Olympic Elite Exam Lift Table.




Figure 2

Item	Description	Item	Description
A	Total-Lock Casters	E	Foot switch Guard
B	Power Strip	F	Telescoping enclosure panels
C	Power Cord	G	Corian® table top
D	Foot switch (Identical switch on opposite side)	H	Scale (optional)

Introduction

The Olympic Elite Exam Lift Table is an adjustable-height examination table designed to provide better patient presentation, save time and improve operator comfort. It is suitable for examining patients up to 250 pounds anywhere within the clinic. Height is fully adjustable from 16" to 39" so heavy patients can be easily loaded from the floor and raised to the optimal height for examination. The Corian® tabletop provides a warm, patient-friendly surface that is easy to clean and always looks new. The built-in scale (optional) provides live weights at any time with the push of a single button.

Safety

Please read and follow the instructions in this manual, including all safety instructions and precautions, carefully and completely.

The user should be aware that the Exam Lift Table is a powerful machine with moving parts capable of causing injury if proper precautions are not taken.

Caution: The lift is raised and lowered by means of an electrical actuator, and care must be taken to assure that hands, feet and other body parts are not in the path of the telescoping panels as the unit is lowered.

Caution: Do not use the Exam Lift Table for patients weighing more than 250 pounds.

Unpacking and Setup

Warning: The Olympic Elite Exam Lift Table is heavy and requires at least 2 people to unpack and move in order to avoid injury to staff or damage to the product.

Note: If you suspect concealed shipping damage, please save all packaging and contact Olympic Veterinary immediately at 1-800-789-6867. Please inspect product carefully after putting it into use and notify Olympic Veterinary within 10 days regarding any initial quality issues.

1. Move the pallet to a level location as close as possible to the intended installation site.
2. Remove the shrink wrap and the cardboard lid, then completely peel off the cardboard shroud.
3. The table is bolted to the floor of the pallet with 4 L-shaped brackets **A**, each secured by 3 hex head screws. Remove all 12 screws using the included hex head wrench.
4. Remove the 2x4 retaining board **B** located in front of the exam lift table from the floor of the pallet by sliding it out of the slotted enclosure.
5. Position the two included wood ramps **C** in front of the exam lift table so that the protruding pin **D** at the end of each ramp drops into the matching hole in the floor of the pallet. Test ramp to assure it is securely attached.
6. Unlock all 4 casters **E** and then slowly lower the table down the ramp until it is on level ground. It can now be rolled to its final destination.
7. Remove the cardboard box stored in the lower deck of the pallet **F**. This contains the Corian® top, which can now be placed over the metal tabletop of the exam lift table.




Figure 1

Assembly

The only assembly required is to lower the Corian® top onto the integrated metal top. The top will fit snugly and does not require fastening. The scale display (optional) is pre-calibrated at the factory. Simply install the AA batteries (provided) and it is ready to use.

Operation

Loading/Unloading Patients

Caution: Be sure the Exam Lift Table is located on a level surface and that all 4 casters are locked before loading or unloading patients.

Caution: Do not use the Exam Lift Table for patients weighing more than 250 pounds.

1. Lower the table top to the desired height for loading a patient by pressing and holding the up or down indicator on either of the two foot switches.
2. Place the patient on the center of the table and hold patient while table is lowered or raised.

Height Adjustment

The height of the Olympic Elite Exam Lift Table is infinitely adjustable from 16" to 39" and may be adjusted at any time during an examination to provide maximum comfort to the operator or to provide better patient visualization.

Two foot switches are provided for convenient access from either side of the table. Care should be taken to keep the toe of operator's shoe out of the path of the descending telescoping panels. A built-in toe guard encourages safe use.

Caution: Care must be taken to assure that hands, feet and other objects are not in the path of the telescoping panels as the unit is lowered.

Scale

The Olympic Elite Exam Lift Table (model 51202) includes a precision 250 lb. capacity scale. The Olympic Scale comes pre-calibrated, so there is no setup procedure. Just install 6 AA batteries (provided) and it is ready to use.

Operating Procedure:

1. **Press WEIGH** – Turn the scale on by pressing the **WEIGH** button. With nothing on the table top, it will display zero (0.0). If there is a patient on the table top, the scale will immediately begin displaying a series of weights as it seeks to find a stable average weight for the patient.

Note: Should the scale display something other than Zero (0.0) with no patient on the table, press and hold the **WEIGH** button for 3 seconds to re-zero the scale. This would be necessary, for example, if a rubber mat or some other external object has been added to the table top for this weighing.

2. **Wait 3 seconds for averaged weigh to lock on display** – Provided the patient is not moving excessively, the scale will auto-lock on an averaged weight within 3 seconds. A flashing arrow will appear in the display window while the scale is averaging. A lock icon will appear in the window when a final weight has been determined. If necessary, the patient can be reweighed by pressing the **WEIGH** button again.




Figure 3

Note: The operator can over-ride the averaging feature at any time by momentarily pressing the **HOLD** button. This will lock on the instantaneous weight then displayed. This feature might be useful, for example, on a very active animal where the scale cannot find a stable average weight.

- 1. Press LB/KG Units button.** Allows the user to toggle between pounds and kilograms at any time. Lb or Kg will appear in the display window.
- 2. Scale automatically turns off.** The scale will automatically turn off 60 seconds after the last button is pushed or after last weight is locked on the display. The last displayed weight can be retrieved at any time by pressing and holding the dual-purpose **HOLD/RECALL** button for 3 seconds.

Note: To TARE out the weight of any foreign object that will remain on the scale during weighing, place this object on the table top and press and hold the **WEIGH** button for 3 seconds. This will provide a new zero (0.0) that has taken out the weight of the foreign object. The patient is now weighed as described above. Care should be taken to remove the TARE at the conclusion of this weighing, by removing all foreign objects and pressing and holding the **WEIGH** button again for 3 seconds until 0.0 is displayed.

Power

The Olympic Elite Exam Lift Table operates off 120V/60Hz line voltage. A 4 port power strip with a 10' power cord is mounted on the end of the unit under the base plate. The lift motor is permanently plugged into one of these ports. The other 3 ports are available as needed.

Moving the Elite Exam Lift Table

If at any time you wish to move the table, disconnect the Power Cord and unlock all four caster brakes. To unlock a caster, step on the brake release lever.

Cleaning

The Corian® top and all painted surfaces – including the base, telescoping column, and foot switches can be cleaned with soap and water and ordinary disinfectant.

Routine Maintenance

The Exam Lift Table requires no routine maintenance.

For any questions about the use or maintenance of the Table, please call 1-800-789-6768.

Specifications

Olympic Elite Exam Lift Table – Model 51201, 51202

Capacity – 250 pounds with centered load

Table size – 21"W x 54"L

Overall maximum dimensions – 21"W x 62"L x 29"H

Power – 120V/60Hz

Height range – 16" to 39"

Model 51201– Olympic Elite Exam Lift Table with scale.

Model 51202 – Olympic Elite Exam Lift Table without scale.

3-Year Unconditional Warranty

The Olympic Elite Exam Lift Table (model 51201, 51202) is unconditionally warranted for a period of three (3) years from the date of invoice.

Olympic Veterinary will elect to repair or replace a defective component without charge within the warranty period. The charge for shipping the defective component is the responsibility of the buyer and must be accompanied by an RMA number. The shipping charge to return the component to the buyer is the responsibility of Olympic Veterinary, Inc. In some cases, at Olympic Veterinary's sole discretion, it may be possible to provide the user with a repair kit, authorizing the user's technical staff to complete the repair. In such cases, the labor cost shall be the sole responsibility of the user.

This warranty is not valid for damage resulting from negligence, accident, unreasonable use, abuse or misuse, exceeding weight limits or altering the product without prior Olympic Veterinary authorization.

Olympic Veterinary expressly disclaims and excludes any liability for consequential, incidental, indirect or punitive damages or financial loss to people or property resulting from any breach of warranty or from the improper operation of this equipment.

This limited warranty is provided to the original purchaser and is not transferable. This constitutes the complete and final agreement for Olympic Veterinary's Elite Exam Lift Table, model 51201, 51202.

Contact Us

Olympic Veterinary Corporation, Customer Service

9311 SE 36th Street, Unit 220
Mercer Island, WA 98040

Phone: (800) 789-6867

Fax: (800) 905-0839

Email: info@OlyVet.com

NOTE: THIS IS NOT THE STREET ADDRESS TO BE USED FOR RMA RETURNS.

Detailed shipping instructions will be provided by customer service when a RMA is issued.